

# Know before use!

## 1. Ramp cleaning

Before using the plasma lamp and rolling ball, lightly wipe them with a wet tissue/dry cloth before use.



## 2. Battery charging

If battery discharge continues, it will affect the lifespan of the device, so keep it charged.



## 3. Power button

The button is located only in the center, lightly press the area below the LED.



## 4. charging cap

This product uses high voltage, so wear a charging cap to avoid touching it with your hands.



# Rolling container combination

## 01. Injecting cosmetics 2. Separating the ring

After opening the rolling device, insert cosmetics only up to the inner thread line (25ml).

※ Avoid getting cosmetics on the screw line.

If you slightly turn the 'container coupling ring' clockwise, the ring will separate with a 'click' sound



## 3. Combine

After inserting the rolling container, fit the protrusions on the container coupling ring into the grooves of the device.

## 4. Fixing

Secure it with counterclockwise force until you hear a 'click' sound.

Shake in the downward direction 1-2 times before use.

